


# SoRTEV™ RNA Enrichment Kit

User Manual – Version 2020/12 Rev.05

For Research Use Only



24 reactions



EXO-SOR-LV



**EXOSOMICS SpA**

Address: Strada del Petriccio e Belriguardo 35, 53100 Siena, Italy

Tel: +39 0577 381308/ Fax +39 0577 381452

Business Registry: SI-137798 VAT#: IT01317090528

[info@exosomics.eu](mailto:info@exosomics.eu)

[www.exosomics.eu](http://www.exosomics.eu)

## INDEX

Intended Use .....	3
Product Description .....	3
Materials Provided/Required .....	3
Kit Components .....	3
Materials Required but Not Provided .....	4
Method Description and Procedure .....	4
Method Description.....	4
Procedure .....	4
Trouble-shooting .....	8
Warnings and Precautions.....	8
Safety Information.....	8
General Precautions .....	9
Limitations .....	9
Symbols.....	9
Contact Information .....	9

### Intended Use

**SoRTEV™ RNA Enrichment kit** is a pre-analytical kit to isolate RNA contained in tumor enriched exosomes from plasma.

The **SoRTEV™ RNA Enrichment Kit (EXO-SOR-LV)** kit is optimized for input volumes ranging from 0.5 ml up to 2 ml of plasma.

The kit is designed to be used with any downstream application employing enzymatic amplification or enzymatic modifications of RNA isolated with SoRTEV, followed by signal detection or amplification (PCR, real time PCR and digital PCR).

Exosomics **SoRTEV™ RNA Enrichment kit** does not provide a diagnostic result. It is responsibility of the user to use and validate the kit in conjunction with any downstream assays.

### Product Description

**SoRTEV™ RNA Enrichment kit** is ready-to-use and it is meant for running 24 tests. Kit allows the isolation of nucleic acids from tumor enriched exosomes, from a minimum of 500 µl of plasma following two subsequent working steps:

1. Tumor-exosome isolation from biofluid.
2. RNA purification.

The purification is based on Exosomics proprietary affinity method and does not require any special equipment, such as ultracentrifugation or chromatography, with a turnaround time of a minimum of 3 hours.

### Materials Provided/Required

Kit components, meant to run a total of 24 reactions, and their storage conditions are listed below.

### Kit Components

Product Code			EXO-SOR-LV (0.5-2ml)
Components	Name	Description	Quantity (number/volume)
<b>Isolation Agent</b>	IA-SOR	Reagent for isolation	1 vial (100 µl)
<b>10X Isolation Buffer</b>	IB-SOR	Diluent for isolation	1 bottle (15 ml)
<b>10X Bead Wash Buffer</b>	BW-SOR	Buffer for bead washing	1 bottle (15 ml)
<b>Isolation Tubes</b>	EXO-IsoT-2ml	Tubes for EV isolation	48 tubes (2 ml)
<b>Lysis Buffer</b>	LB-SOR	Buffer for vesicle lysis	1 bottle (20 ml)
<b>Washing Buffer</b>	WB-SOR	Buffer for column washing	1 bottle (30 ml)
<b>Elution Buffer</b>	EB-SOR	Buffer for RNA elution	1 bottle (4 ml)
<b>RNA purification columns</b>	RC-SOR	Columns for RNA purification	24 columns
<b>Elution Tubes</b>	EXO-CoIT-1.5ml	Tubes for pure RNA collection	24 collection tubes (1.5 ml)

Store all the supplied reagents according to the instructions on the respective packages. In particular:

- Upon arrival store IA-SOR, BW-SOR bottles at +2/+8 °C;
- Upon preparation of WB-SOR reagent, store buffer at +2/+8 °C;
- If properly stored, all the reagents provided with the kit are stable until the expiration date printed on the labels.

### **Materials Required but Not Provided**

- 96-100% Ethanol (Example: RNAase-free Ethanol Sigma Aldrich, Cat. Num. 51976)
- Chloroform (Example: Sigma Aldrich Cat. Num. C2432)
- RNAse-free 2 ml tubes for Molecular Biology (24 tubes, Example: Starlab Cat. Num. S1620-2700)
- Disposable Gloves
- Single-use and/or pipettes with disposable tips
- Pipettes for reagent preparation
- Ultrapure water for dilution of (10X) **IB-SOR** and (10X) **BW-SOR**
- Ice for EV lysis during RNA purification
- Benchtop centrifuge with rotor for 2 ml reaction tubes (kit validated on Eppendorf 5415R)
- Rotator or tilting shaker
- Vortex

### **Method Description and Procedure**

#### **Method Description**

**Method:** SoRTEV™ RNA Enrichment Kit method is based on Exosomics' proprietary immuno-affinity method that enriches for tumor derived exosomes from which EV-RNA is then extracted. SoRTEV™ RNA affinity method yields the highest level of enrichment of tumor genetic material.

**Sample volume:** SoRTEV™ RNA kit has been optimized for sample volumes ranging from 0.5 ml to 2 ml of plasma. Follow **protocol A** up to 1 ml of plasma. For volume >1 ml of plasma, the best performance is obtained by splitting plasma into two 2 ml vials (EXO-IsoT-2ml), as described in **protocol B**.

**Sample collection:** Exosomics S.p.A. suggests the following recommendations,

- tubes for blood collection with anticoagulant, K2 EDTA, or CTAD or citrate;
- store the collection tubes at room temperature (20-25°C);
- process the samples within 30 minutes or in a shorter time possible after collection;
- centrifuge the blood samples at 1500 g 10 minutes at room temperature (20-25°C);
- recover the plasma on the top, being careful to not contaminate the samples with red blood cells;
- aliquoting is recommended since freeze-and-thaw cycles reduce the quality of the sample;
- store the plasma collected freeze (preferred temperature -80°C) or proceed with the SoRTEV™ RNA Enrichment Kit protocol.

#### **Procedure**

SoRTEV™ RNA allows the isolation of tumor originated exosomes from plasma, and extraction of their RNA through a four-step procedure:

1. Reagent preparation
2. Plasma preparation
3. EV isolation from plasma
4. EV-RNA purification

## 1 Reagent preparation:

Equilibrate all reagents to room temperature before the use.

- 1.1 **1X Isolation Buffer (1X-IB):** dilute 10X Isolation Buffer (**IB-SOR**) in fresh Ultrapure water to a final 1X concentration (i.e. 1 ml of IB-SOR and 9 ml of ultrapure water) and label the vial as “1X-IB”.
- 1.2 **1X Bead Wash Buffer (1X-BW):** dilute 10X Bead Wash Buffer (**BW-SOR**) in fresh Ultrapure water to a final 1X concentration (i.e. 1 ml of BW-SOR and 9 ml of ultrapure water) and label the vial as “1X-BW”.
- 1.3 **Washing Buffer (WB):** add 20.9 ml of pure Ethanol (96-100%) in **WB-SOR** bottle. Mix well by inverting 6-8 times.

## 2 Plasma preparation:

### Protocol A. Plasma volume ≤ 1 ml

- 2.1a Pre-clear the plasma sample by centrifuging at 1200 g for 20 min at 10 °C to eliminate red blood cells and cellular debris.
- 2.2a Transfer pre-cleared plasma in the appropriate tube (EXO-IsoT-2ml, included in the kit) and discard the pellet and debris.
- 2.3a Dilute pre-cleared plasma in 1:1 v/v with 1X-IB (refer to **Table 1** for details; e.g. if processing 0.5 ml of plasma, add 0.5 ml of 1X-IB).
- 2.4a Proceed to EV isolation from plasma with **protocol A**.

### Protocol B. plasma volume > 1 ml

- 2.1b Pre-clear the plasma sample by centrifuging at 1200 g for 20 min at 10 °C to eliminate red blood cells and cellular debris.
- 2.2b Transfer pre-cleared plasma in the appropriate tubes (EXO-IsoT-2ml, included in the kit) splitting the available volume in two tubes equally (e.g. if processing 2 ml of plasma, transfer 1 ml plasma in the first tube and 1 ml in the second one). Discard the pellet and debris.
- 2.3b Dilute pre-cleared plasma in 1:1 v/v with 1X-IB (refer to **Table 2**; e.g. if processing 2 ml of plasma, add 1 ml of 1X-IB in the first tube and 1 ml in the second one).
- 2.4b Proceed to EV isolation from plasma with **protocol B**.

Plasma Volume (ml)	1X-IB to be added (ml)	Total sample volume (ml)	Quantity of EXO-IsoT-2ml tubes	IA-SOR to be added (µl)
0.5	0.5	1	1	2
0.8	0.8	1.6	1	2
1	1	2	1	2

**Table 1: Plasma volume ≤ 1 ml.** The table summarizes how to perform plasma dilution, how many isolation tubes (EXO-IsoT-2ml) to use, and the volume of isolation agent (IA-SOR) to be added per reaction.

Plasma Volume (ml)	1X-IB to be added (ml)	Total sample volume (ml)	Quantity of EXO-IsoT-2ml tubes	IA-SOR to be added (µl)
1.5 (0.75 + 0.75)	1.5 (0.75 + 0.75)	3	2	2 + 2
2 (1 + 1)	2 (1 + 1)	4	2	2 + 2

**Table 2: Plasma volume > 1 ml.** The table summarizes how to perform plasma dilution, how many isolation tubes (EXO-IsoT-2ml) to use, and the volume of isolation agent (IA-SOR) to be added per reaction.

### 3 EV isolation from plasma

#### Protocol A. Plasma volume ≤ 1 ml

- 3.1a Add 2.0 µl of antibody-coated beads reagent (IA-SOR) to the pre-cleared diluted sample.
- 3.2a Mix well by inverting the tube/s.
- 3.3a Incubation time is 2 hours at RT under rotation.
- ⓘ If a rotator is not available, a tilting shaker can be employed for sample mixing during incubation.
- 3.4a Centrifuge 10 minutes at 9300 g at RT.
- 3.5a Discard the supernatant and resuspend the pellet by gently adding 1 ml of 1X-BW.
- 3.6a Spin the sample at 9300 g for 10 min at RT.
- 3.7a Discard the supernatant again and resuspend the pellet by gently adding 1 ml of 1X-BW.
- 3.8a Spin the sample at 9300 g for 10 min at RT.
- 3.9a Discard the supernatant.
- 3.10a Proceed to EV Lysis (4.1).

#### Protocol B. Plasma volume > 1 ml:

- 3.1b Add 2.0 µl of antibody-coated beads reagent (IA-SOR) into each (2) isolation tube (EXO-IsoT-2ml) and follow the procedure using two tubes. Refer to **Table 2** for details.
- 3.2b Mix well by inverting the tubes.
- 3.3b Incubation time is 2 hours at RT under rotation.
- ⓘ If a rotator is not available, a tilting shaker can be employed for sample mixing during incubation.
- 3.4b Centrifuge 10 minutes at 9300 g at RT.
- 3.5b Discard the supernatant and resuspend the pellet by gently adding 1 ml of 1X-BW.
- 3.6b Spin the sample at 9300 g for 10 min at RT.
- 3.7b Discard the supernatant again.
- 3.8b Pool the two pellets of the same sample in one vial: resuspend pellet 1 by gently adding 1 ml of 1X-BW and transfer the resuspended pellet into the vial containing pellet 2, also resuspending the second pellet. You have pooled the two pellets in the same vial. Discard the tube that contained pellet 1.
- 3.9b Spin the sample at 9300 g for 10 min at RT.
- 3.10b Discard the supernatant.
- 3.11b Proceed to EV Lysis (4.1).

## 4 RNA purification:

### 4.1 EV Lysis:

4.1.1 Add 700 µl of lysis buffer (LB-SOR) and vortex 30 seconds.

4.1.2 Incubate 5 minutes at RT.

**i** At this stage it is possible to freeze the sample at -80°C.

4.1.3 Add 140 µl of pure chloroform (not provided with the kit).

4.1.4 Shake the tube for 30 seconds.

4.1.5 Incubate 10 minutes at RT.

4.1.6 Incubate 1 minute in ice and centrifuge at 12000 g at 4 °C for 10 minutes.

4.1.7 Transfer the top phase in a fresh tube (RNase-free 2 ml tube, not provided with the kit).

4.1.8 Add ethanol (96-100%) to the recovered phase in a 2:1 v/v ratio (i.e. add 640 µl of ethanol to 320 µl of recovered phase). Mix gently inverting 4-5 times.

### 4.2 RNA purification:

4.2.1 Transfer 600 µl of the mixture into a spin column (RC-SOR).

4.2.2 Spin at 14000g for 30 seconds at RT.

4.2.3 Discard the flow through.

4.2.4 Repeat steps 4.2.1-4.2.3 with the remainder.

4.2.5 Add 400 µl of RNA Washing Buffer (WB-SOR) in the spin column (RC-SOR).

4.2.6 Spin at 14000 g for 30 seconds at RT.

4.2.7 Discard the flow through.

4.2.8 Repeat steps **4.2.5 - 4.2.7** two more times for a total of three washes.

4.2.9 Spin 2 additional minutes at 14000 g at RT to eliminate ethanol residues from the column.

4.2.10 Remove the tube and transfer the spin column into an elution tube (EXO-CoIT-1.5 ml).

4.2.11 Elute the column with 15 µl of elution buffer (EB-SOR).

4.2.12 Incubate 5 minutes at RT.

4.2.13 Spin 2 minutes at 200 g at RT. Spin 1 minute at 14000 g at RT, keep the flow through.

4.2.14 Eluted RNA is now ready for downstream analysis or for storage at -80 °C.

**i** For low abundant targets, we advise to proceed immediately to downstream analysis, to avoid RNA degradation during freeze-thawing cycles.

## Trouble-shooting

This table may solve some technical problems that could arise during **SoRTEV™ RNA Enrichment Kit** protocol execution. For more information, please contact us at [support@exosomics.eu](mailto:support@exosomics.eu).

Technical Problems	Potential Causes	Suggestions and comments
Low RNA recovery	Poor plasma quality due to delayed blood processing	Repeat blood processing according to Sample processing in “Method Description and Procedure” section.
	Plasma samples are frozen and thawed multiple times	Always use fresh samples or samples thawed once.
	Prolonged sample storage at room temperature	Do not keep the samples at RT for prolonged time.
	Wash Buffer (WB-SOR) prepared incorrectly	Check that these buffers were diluted in the correct volume of 96-100% ethanol.
	The eluate volume is lower than the applied volume	Expect to recover an eluate volume with 2-3 µl less than the applied volume due to retention of the silica membrane.
RNA not suitable for enzymatic reaction	Presence of ethanol traces in eluate	Make sure to remove all ethanol residuals from the column (RC-SOR) before eluting the sample.
	RNA degradation	Avoid RNA freeze-thawing cycles and keep it on ice while working. For long term storage, keep it at -80°C.
Presence of DNA traces in RNA eluate	DNA carry over during RNA purification	DNase treatment can be performed (Example: Sigma-Aldrich, Cat. Num. AMPD1-1KT) on beads/EVs pellet (after <b>step 3.9a or 3.10b</b> ). If using Sigma-Aldrich Amplification Grade DNase I, incubate sample with DNase I and 1X Reaction Buffer for 15 min at RT. The reaction does not require inactivation, only addition of lysis buffer is required.

## Warnings and Precautions

**SoRTEV™ RNA** Enrichment kit does not provide a diagnostic result. It is the sole responsibility of the user to use and validate the kit in conjunction with a downstream assay. To minimize irregularities in results, suitable controls for downstream applications should be used.

All products sold by Exosomics are subject to extensive quality control procedures and are warranted to perform as described and used correctly. Any problems should be reported immediately. The kit contents are for laboratory use only, and they must be stored in the laboratory and must not be used for purposes other than intended.

## Safety Information

When working with chemicals, always wear a suitable lab coat, disposable gloves, and protective goggles. Lysis Buffer (**LB-SOR**) contains phenol and guanidinium salts, and should be handled with care. Phenol should be handled under chemical hood. Guanidinium salts form highly reactive compounds when combined with bleach, thus care must be taken to properly dispose of any of these solutions.

For more details, please refer to the **SoRTEV™ RNA Enrichment Kit** Safety Data Sheet.



## General Precautions

Always operate in accordance with good laboratory guidelines and the instructions included in this Handbook.

Handle and dispose of waste materials and reagents from the use of this kit in accordance with current regulations on staff safety and environmental protection.

## Limitations

Results from the product must be interpreted within the context of all relevant clinical and laboratory findings and are not to be used alone for diagnosis.

The product is to be used only by professional users such as technicians, physicians and biologists experienced and trained in molecular biological techniques, including RNA isolation.

Dilution of the reagents, other than as described in this handbook, is not recommended and will result in a loss of performance.

Attention should be paid to expiration dates and storage conditions printed on the box and labels of all components. Do not use expired or incorrectly stored components.


## Symbols



Contains reagents sufficient for n tests



Catalogue number



Batch or lot code



Use by



Temperature limit



Operating instructions



Manufacturer

## Contact Information

**Exosomics S.p.A.**

**Address** Strada del Petriccio e Belriguardo 35, 53100 Siena, Italy

**Business Registry:** SI-137798 **VAT#:** IT01317090528

**Administration:**

Tel: +39 0577 381308

For info: [administration@exosomics.eu](mailto:administration@exosomics.eu)

**Technical Support:**

For info: [support@exosomics.eu](mailto:support@exosomics.eu)

**Sales & Marketing:**

Tel: +39 0577 381309